




**Kalika Manavgyan Secondary School  
Butwal-10-Rupendehi**



**Project Work  
Subject: Chemistry**

**Detection of Presence of Acetic Acid In Vinegar**

# Acknowledgement

I hereby declare that the project work “ **Detection of Acetic Acid in Vinegar**” submitted to Department of Chemistry **Kalika Manavgyan Secondary School** in the form of hard copy of project work which has done under the supervision of Chemistry teachers **Kabir Oli** and **Amrit Khanal** and is submitted for the partial fulfillment of the requirements for the secondary level education of National Examination Board Chemistry Grade 12.

## Group Members

1. Shivpal Chaudhary
2. Sushant Bhattra
3. Rajan Paudel
4. Rupesh Baral
5. Shiva Neupane
6. Samrat Basual

# Index

★ Object.....	06
★ Background of the study .....	07
★ Materials and Chemicals .....	08
★ Theory.....	09
★ Experimental Procedure.....	10
★	
★ Observation.....	11
★ Result Conclusion and Precautions .....	11
★ Bibliography.....	11

**TO DETECT  
THE PRESENCE  
OF ACETIC  
ACID IN A  
SAMPLE OF  
VINEGAR**

## BACKGROUND OF THE STUDY

Vinegar is a solution made from fermentation of ethanol. There are many types of vinegar each starting from different sources. The amount of acetic acid( $\text{CH}_3\text{COOH}$ ) in vinegar can vary, typically between 4% to 6%.

In this project, we determine the amount of acetic acid in different vinegar using titration. Titration is a way to measure the unknown amount of a chemical(titrant) in a solution with a known concentration(titrant) and the end point is monitored.

To monitor the acidity of a vinegar solution we can add enough hydroxyl ions to balance out the added hydrogen ions from the acid. Hydroxyl ions react with hydrogen ions to give water.

For titrand we use a dilute solution of sodium hydroxide. Sodium hydroxide is a strong base and it disassociates into sodium and hydroxyl ions in water.

In this experiment we will use phenolphthalein as indicator.

Phenolphthalein is colourless when the solution is acidic or neutral

When the solution becomes slightly basic, phenolphthalein turns pinkish.

## **APPARATUS REQUIRED**

- 1.Small funnel
2. Conical flask
3. 10ml Burette
4. Stand
5. Pipette
6. Burette Clamp

## **CHEMICALS REQUIRED**

1. Vinegar
2. Water
3. Sodium Hydroxide Solution
4. Phenolphthalein

## THEORY


Required amount of sodium hydroxide can be calculated by formula:

$$W = \frac{\text{Molarity} \times \text{Molar mass} \times \text{Volume (cm)}^3}{1000}$$

Molar mass of NaOH = 40 g/mol

$$W = \frac{0.5 \times 40 \times 500}{1000} = 10 \text{ g}$$

The acetic acid content of vinegar may be determined by Titrating vinegar sample with known molar concentration(molarity).



At the endpoint of titration, according to law of stociometry

$$N_1 \times V_1 = N_2 V_2$$

where,

N1 = Concentration of acid

N2 = Concentration of base

V1 = Volume of acid

V2 = Volume of base

Finally, we can use the normality found to calculate the concentration of acetic acid in gram per liter as:

$$\text{gram/liter} = \text{Normality}(N_1) \times \text{Equivalent weight}$$

=

$$= N_1 \times 60 \text{ g/l}$$

## **PROCEDURE**

1. We poured vinegar in a conical flask.
2. We added water to dissolve the vinegar so that the volume of solution becomes 20ml.
3. We added 2 drops of phenolphthalein solution.
4. We used burette clamp to attach the burette to the ring stand.
5. We used funnel to fill the burette with 0.5N solution of Sodium hydroxide.
6. We noted the starting level of sodium hydroxide solution in burette.
7. We put vinegar solution to be titrated below the mouth of burette.
8. We slowly dripped the solution of sodium hydroxide into the vinegar solution
9. When the vinegar turned pink and the color remained even after mixing, we quickly noted the level of sodium hydroxide in the burette after closing the tap.
10. We noted the remaining level of sodium hydroxide by reading from bottom meniscus.



=

11. We subtracted the initial level from remaining level to figure the volume of NaOH used.

12. We repeated the process for 3 times.

## OBSERVATION

The following table shows the outcomes of the titration.

S No.	Volume in Vinegar Solution (in ml)	Burette Reading		Volume of NaOH solution used (in ml)
		Initial (in ml)	Final (in ml)	
1	20	0	15.2	15.2
2	20	15.2	30.4	15.2
3	20	30.4	45.6	15.2

We take 15.2 ml for volume of NaOH during calculation as it is concordant.

## **CALCULATION**

We know that

$$\begin{aligned} N_1 \times V_1 &= N_2 \times V_2 \\ \text{or } N_1 \times 20 &= 0.5 \times 15.2 \\ \text{or } N_1 &= 0.38 \text{ mol/L} \end{aligned}$$

Now we calculate,

$$\begin{aligned} \text{concentration of acetic acid} &= N_1 \times 60 \\ &= 22.8 \text{ g/litre} \end{aligned}$$

## **RESULT**

So, the strength of acetic acid in a sample of vinegar is found to be 22.8 gram per litre.

## **CONCLUSION**

Concentration of acetic acid in vinegar can be calculated by titration in lab.

## **PRECAUTIONS**

1. Wear protective glasses and clothes during experiment.
2. We should not use more than 2-3 drops of indicator as it can alter the volume of solution.
3. NaOH should be dropped slowly drop by drop.
4. We should mix the solution properly
5. We should use clean water for solution

## BIBLIOGRAPHY

[www.wikipedia.org](http://www.wikipedia.org)  
[www.slideshare.com](http://www.slideshare.com)  
[www.byjus.com](http://www.byjus.com)

